

**Indian Council of Agricultural Research
National Institute of High Security Animal Diseases
Anand Nagar, Bhopal – 462021, M.P.**

Electron Microscopy Unit

Proforma for sample submission by external users

For Office Use	
Requisition Number	
Accession Numbers	

Date of submission :

Name, designation and address of the investigator :

Telephone: Email :

Name, designation and address of the user :

Telephone: Email :

Name of the organization/institute :

Category of organization : I / II*

Name of the project or title of dissertation or thesis :

Nature of samples :

Details of fixative or diluents used :

Number of specimen to be processed :

Details of Demand Draft : No.....Date.....

Bank:

Billing address :

Signature of the investigator
Date and official seal

Note:

- Tissue samples should have been fixed soon after death/sacrifice. Tissues should be cut in to 1mm x 1mm size and fixed in suitable fixative at 4°C for 6 to 12 hours depending on the tissue type. The fixed specimen should be transported in 0.1M phosphate buffer (pH7.4).
- A Maximum of five samples will be accepted at a time.
- * Category I - Private organizations.
Category II - Government/national organizations/educational institutions and universities.